

PHY 517 / AST 443:

Observational Techniques in Astronomy


Lecture 4:

Photometry

Spectroscopy

Please be gentle with our equipment!!

Example: do not force a plug into a connector - USB-A and -B are **not** symmetric. Check with the TAs when in doubt.



For the time being, the STL-100IE camera is out of commission. Plan on using the ST-402ME camera instead.

Shut-Down Procedure

Make sure to follow the entire shut-down procedure in the manual(s).

I found the dome power still on in the daytime; there are open wires in the dome, so this is a safety hazard.

Last Time: CCDs

You want to take a sky image with our telescope + CCD camera. Which calibration images do need to take in order to analyze the science image? How do you take them? What do you do with them?

- **Dark frames:** same exposure time, same temperate -> measures bias + dark current
- **Flat-fields:** images of a uniform surface (e.g. dome), taken with the camera attached to the telescope in the same set-up as for the science image (*also needs to be dark-corrected*)

$$\text{calibrated image} = \frac{\text{science image} - \text{dark}}{\text{flatfield}}$$


Calibrated image

Photometry: measure the brightness (also colors) of sky objects

How would you measure the flux of this star?

1. Sum up counts within aperture
2. Measure sky background (per pixel) from empty region of sky
3. Subtract estimated sky counts
4. Repeat for a reference star of known magnitude

$$m_1 - m_2 = -2.5 \log \left(\frac{F_1}{F_2} \right)$$




Spectroscopy

Motivation

photometry only measures integrated flux (over filter band)

gives some
information
about the object
properties, but
often not
enough

e.g.: finding the most
distant galaxies




Narrow-band imaging

can determine spectrum
of object with images in
many narrow-band filters

advantage: can determine
spectra of all objects in
the same FOV

disadvantage: have to take
a lot of images!




The PAU Survey

Spectroscopy

add a dispersing element
to split up the light from
an object: measure the
spectrum directly

e.g. a prism:




Olympus

Prism Spectroscopy

only few
astronomical
spectrographs
use prisms

- low dispersion (resolution)
- dispersion varies with wavelength



“low dispersion prism” for IMACS spectrograph on Magellan 6-m telescope; uses 3 prisms


Diffraction gratings

make use of wave
properties of light:
interference

grating: many parallel
lines ($\sim 500/\text{mm}$)

similar to single-slit and
double-slit experiments

position of n th order:




D. Pogosian

if $b-c = \lambda$: maximum at θ
and $d-e = 2\lambda$, etc.

$$n\lambda = d \sin \theta$$


Diffraction gratings



Diffraction gratings

can be transmission gratings or reflection gratings

most astronomical spectrographs use reflection gratings





Dhillon 2012

blaze wavelength: wavelength for direction of reflection coincides with desired spectral order
→ maximal efficiency

Typical spectrograph

entrance: usually a slit, similar to seeing size

collimator: converts a diverging beam to a parallel beam



A Schematic Diagram of a Slit Spectrograph

Spectral Resolution

defined by smallest wavelength difference $\Delta\lambda$ that can be distinguished at wavelength λ

$$R = \frac{\lambda}{\Delta\lambda}$$

determined by:

- grating (line density)
- width of entrance slit
- seeing

resolution: R or $\Delta\lambda$

dispersion: length $\Delta\lambda'$ of spectrum over single pixel, [$\text{\AA}/\text{px}$]

to properly sample the spectrum:

$$\Delta\lambda \sim 2 - 3 \Delta\lambda'$$




Spectral Resolution

$R < 1000$	low-resolution	e.g. our “low-resolution” spectrograph
$1000 < R < 10,000$	medium-resolution	e.g. our “high-resolution” spectrograph
$R > 10,000$	high-resolution	Echelle spectrographs

Long-slit observations


Planetary Nebula Spectroscopy : NGC 6720 [Ring Nebula] & NGC 6543 [Cat's Eye Nebula]

Jim Ferreira, Livermore CA



Long-slit observations

D. Steeghs




- long axis of CCD used to sample spectrum
- spatial information along slit still available: two objects, lots of sky
- sky background has a lot of emission lines!

Spectroscopic Calibration

- dark frames!
- flat field: use bright continuum source
 - small-scale pixel sensitivity variation
 - variations in slit width
- wavelength calibration: which position on the CCD corresponds to which wavelength?
 - use “arc” lamps with discrete emission lines
 - can also use sky emission lines
- flux calibration:
 - “spectrophotometric” standard stars: stars with known spectral shapes, smooth continua

Spectroscopic Calibration

wavelength calibration: map pixel position to emission lines




C. Buil /
astrosurf

Spectroscopic Calibration

flux calibration:
observe
spectrophotometric
standard star

compare observed
spectrum (counts)
to known spectrum




Long-slit spectrographs

most common spectrograph

can only target one (or a few) objects


gives spatial variation

very good estimate of sky background



Echelle spectrographs

- very high resolution long-slit spectrographs
- have additional elements to fit entire spectrum onto CCD
- only for bright objects




Echelle spectrum
of the Sun,
4000-7000Å


NOAO

Multi-object spectrographs

make a mask with multiple slits, one per target




Open Cluster NGC 330 in SMC - VLT UT1 + FORS1 (MOS-mode)



Spectra of Stars in Open Cluster NGC 330 in SMC - VLT UT1 + FORS1 (MOS-mode)

Integral-Field Units


divide image into “spaxels” (spectroscopic pixels)



M. Westmoquette


Lab 2

- Option A) Exoplanet transit
 - photometry of a star that hosts an exoplanet
 - extensive lab manual
 - usually done with STL-1001E ...
- imaging data analysis is easier (conceptually) than spectroscopic data analysis
- involves more programming later on
- for Lab 3, have to then use spectroscopy



Lab 2

- Option B) Spectroscopy of bright stars
 - measure strength of absorption line features as function of temperature
 - lab manual still in draft form
 - spectroscopic data analysis conceptually more involved than imaging data analysis
 - somewhat less programming intensive
 - for Lab 3, have to then do imaging



Lab 2

- Start discussing which Lab 2 option you want to choose
- Should have completed at least the day-time observations of Lab 1 before you decide
- Will look to scheduling these between Feb. 14 and March 11
- Once you have decided, let us know. THEN: select targets according to instructions and request 3 observing nights - first come, first serve